

MolDan Service

REPLACEMENT CT TUBES

MAINTAIN YOUR CT'S TECHNICAL PERFORMANCE

Dunlee replacement CT tubes offer excellent quality and are meticulously engineered to be OEM-equivalent.

To help you maintain your CT's performance, we offer fast delivery of quality CT tubes. We provide the **24/7/365 technical support** you need to replace, install and calibrate your replacement CT tube.

Call +373 22 80 80 22 today or visit www.moldanservice.md to learn more

REDUCE REPLACEMENT COSTS AND SCANNER DOWNTIME

For independent service organizations (ISOs) and in-house teams replacing CT tubes at hospitals and imaging centers, Dunlee's replacement tubes and support services help to reduce replacement costs and scanner downtime.

EXPERIENCE FOR THE FUTURE

Building on 100+ years of experience, Dunlee's people, products, programs and partnerships have shaped how diagnostic images are made. Dunlee remains at the forefront of medical imaging as an international leader in the research, design, and manufacturing of high-performance CT tubes.

TubeType		Ordering No.	Focal Spot	kHU	Standard Warranty
General Electric					
REEVO Series					
	REEVO240G This 8.0 MHU CT tube is designed as a replacement for use in GE Discovery™, GE Lightspeed™ and GE Optima™*	9896-055-96401	0.7 x 0.6 0.9 x 0.9	8,000	12 months, full
DA200 Series					
	DA200P40 This 6.3 MHU CT tube is designed as a replacement for use in GE Optima™*	4598-006-53451	0.9 x 0.7 1.2 x 1.1	6,300	12 months or 6,000 patient exams, full
	DA200ULTRA This 6.3 MHU CT tube is designed as a replacement for use in GE BrightSpeed™, GE Discovery™, GE Lightspeed™ and GE Optima™*	9896-055-88704	0.9 x 0.7 1.2 x 1.1	6,300	12 months, full
DA165 Series		1			
	DA165NP This 3.5 MHU CT tube is designed as a replacement for use in GE HiSpeed, GE BrightSpeed™ and GE Optima™*	9896-055-83901	0.9 x 0.7 1.2 x 1.2	3,500	14,000 slices/90 days fu 140,000 slices or 12 months pro-rated
DA135 Series					_
	DA135 CT/e This 2.0 MHU CT tube is designed as a replacement for use in GE Brivo™, GE HiSpeed and GE CT Dual™*	9896-055-80800	0.7 x 0.6	2,000	8,000 slices/90 days ful 80,000 slices or 12 months pro-rated
Siemens					
S532/Akron Series					
	S532B/Akron B This 5.3 MHU CT tube is designed as a replacement for use in Siemens SOMATOM*	9896-055-95702	0.5 x 0.7 0.6 x 1.3	5,300	100 %: 40,000 ss or 6 months 70 %: 80,000 ss or 8 months 40 %: 150,000 ss or 12 months
	S532Q/Akron Q This 5.3 MHU CT tube is designed as a replacement for use in Siemens SOMATOM*	9896-055-95902	0.5 x 0.7 0.8 x 1.2	5,300	100 %: 30,000 ss or 6 months 70 %: 60,000 ss or 8 months 40 %: 120,000 ss or 12 months

*The products listed may be trademarks of the OEM. For the latest information regarding compatibility for CT tubes and scanners, please refer to our cross-reference guide.

All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Dunlee reserves the right to make changes in specification and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Contact us for further information.

Philips Medical Systems DMC GmbH Tel.: +49 40 34971-1391; www.dunlee.com

Company MolDan Service - Official Distributor in R. Moldova Tel.: +373 22 80 80 22; www.moldanservice.md

CT CROSS REFERENCE

Subject to change without notice. All products listed may be trademarked by the associated OEM.

Brand	CT Scanner Model	Dunlee Replacement Tube	Dunlee Part #	
GE	BrightSpeed™ Edge	DA200ULTRA	9896-055-88704	
GE	BrightSpeed™ Elite	DA200ULTRA	9896-055-88704	
GE	BrightSpeed™ Excel	DA200ULTRA	9896-055-88704	
GE	BrightSpeed Edge™ (Select / Lite)	DA165NP	9896-055-83901	
GE	BrightSpeed Elite™ (Select / Lite)	DA165NP	9896-055-83901	
GE	BrightSpeed Excel™ (Select / Lite)	DA165NP	9896-055-83901	
GE	Brivo™ CT315/325	DA135 CT/e	9896-055-80800	
GE	Brivo™ CT385	DA135 CT/e	9896-055-80800	
GE	Discovery NM/CT570c	REEVO240G	9896-055-96401	
GE	Discovery™ CT590 RT	REEVO240G	9896-055-96401	
GE	Discovery™ NM/CT 670	DA200ULTRA	9896-055-88704	
GE	Discovery™ NM/CT 670 Pro	DA200ULTRA	9896-055-88704	
GE	Discovery™ PET/CT 610	DA200ULTRA	9896-055-88704	
GE	Discovery™ PET/CT 710	DA200ULTRA	9896-055-88704	
GE	Discovery™ RX	DA200ULTRA	9896-055-88704	
GE	Discovery [™] ST (PET/CT)	DA200ULTRA	9896-055-88704	
GE	Discovery [™] STE (PET/CT)	DA200ULTRA	9896-055-88704	
GE	Discovery [™] VCT	REEVO240G	9896-055-96401	
GE	HiSpeed CT/e, CT DualTM	DA135 CT/e	9896-055-80800	
GE	HiSpeed DX/i	DA165NP	9896-055-83901	
GE	HiSpeed FX/i	DA165NP	9896-055-83901	
GE	HiSpeed LX/i	DA165NP	9896-055-83901	





Brand	CT Scanner Model	Dunlee Replacement Tube	Dunlee Part #	
GE	LightSpeed ^{™ 16}	DA200ULTRA	9896-055-88704	
GE	LightSpeed [™] Plus	DA200ULTRA	9896-055-88704	
GE	LightSpeed [™] Pro ³²	REEVO240G	9896-055-96401	
GE	LightSpeed™ QX/i	DA200ULTRA	9896-055-88704	
GE	LightSpeed [™] RT ¹⁶	REEVO240G	9896-055-96401	
GE	LightSpeed [™] Ultra	DA200ULTRA	9896-055-88704	
GE	LightSpeed [™] VCT	REEVO240G	9896-055-96401	
GE	LightSpeed™ VCT Select	REEVO240G	9896-055-96401	
GE	LightSpeed [™] VCT XT	REEVO240G	9896-055-96401	
GE	LightSpeed™ VCT XTe	REEVO240G	9896-055-96401	
GE	LightSpeed™ Xtra	REEVO240G	9896-055-96401	
GE	Optima™ CT540	DA200ULTRA	9896-055-88704	
GE	Optima™ CT580	REEVO240G	9896-055-96401	
GE	Optima™ CT580 RT	REEVO240G	9896-055-96401	
GE	Optima™ CT580 W	REEVO240G	9896-055-96401	
GE	Optima™ CT660	DA200 P40	4598-006-53451	
GE	Optima [™] PET/CT 560	DA200ULTRA	9896-055-88704	
GE	Optima™ PET/CT 560 FX	DA200ULTRA	9896-055-88704	
Siemens	Somatom Sensation 4	S532B/Akron B	9896-055-95702	
Siemens	Somatom Sensation Volume Zoom	S532B/Akron B	9896-055-95702	
Siemens	Somatom Sensation Volume Access	S532B/Akron B	9896-055-95702	
Siemens	Somatom Sensation 10	S532Q/Akron Q	9896-055-95902	
Siemens	Somatom Sensation 16	S532Q/Akron Q	9896-055-95902	
Siemens	Somotom Sensation Cardiac	S532Q/Akron Q	9896-055-95902	

All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Dunlee reserves the right to make changes in specification and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Contact us for further information and to find the tube best suited to your needs.

Philips Medical Systems DMC GmbH Tel.: +49 40 34971-1391; www.dunlee.com

Company MolDan Service - Official Distributor in R. Moldova Tel.: +373 22 80 80 22; www.moldanservice.md